



D TECH DENTAL  
TECHNOLOGIES

## Case Record

Lab Name : Linea Dental Lab, Pondicherry

Case date

12 Nov 2021

Revision No

00

Case No

CR-3D-6/Nov21/C&B/Linea Dental Lab


Revision Date

### 3D Printed Prosthesis

Printed C & B



### Fitment Checking & Reinforcement



Placement In Mouth



### Product used : Accuprint C & B Resin



#### About The case

Patient for full mouth rehabilitation required a 10-unit prosthesis with 5-unit pontic. The 10-unit prosthesis was to be supported on 5 abutments.

#### Outcome & Comment

The prosthesis was printed using D Tech C&B resin. Suitable gap was given between crown and abutment in exocad design software to compensate for shrinkage of 3D printing process. The printed prosthesis was a perfect fit. Some basic polishing techniques were used to achieve a very good level of gloss and aesthetics. Lab was extremely happy with the result.

#### About the Lab

Linea Dental Lab is a high quality lab located in Auroville, Pondicherry. Their master technician is very professional in prosthetic work and specializes in "life like aesthetics"