

 D TECH DENTAL TECHNOLOGIES	Case Record		Case Date	06 Aug 2023
	Doctor Name : Dr.Priyam Aditya, Pune		Case Ref	Aug 23/SG
	Case No	CR-3DR-DB -01	Rev No &Date	02/10.2023

3D Printing - Full Denture Case

Printed Full Denture



Product used : Denture & C&B resin



About The case

Patient came into the clinic with a broken denture and needed an emergency replacement since he had to travel after two days.

Outcome & Doctor's Comment

His existing denture was scanned and the input was taken into Exocad for digital denture design. Base and the teeth were separated using the software feature. Denture base was printed using Accuprint D tech Denture resin and teeth were printed using C&B resin. The two components were assembled and the patient received a fully functional denture within 24 hours. The fitting was exceptionally good and the patient was able to speak clearly and got full functionality. Denture material showed a high ability to be easily polished to a very highly aesthetic level.

About the Doctor

Dr. Priyam Aditya completed his B.D.S. and M.D.S. (Prosthodontics) from BharatiVidyapeeth University, Pune. He is member of Indian Dental Association (IDA), Indian Prosthodontic Society (IPS) and International Team for Implantology (ITI, Switzerland).

