

RESTORE +

High strength posterior Restorative cement

D Tech Product Newsletter



Condensable

Premium quality ,condensable high strength posterior restorative cement for long term restorations with high compressive strength .



Comfort Set

Cement remains in condensable state up to 2 mins from start of mixing, giving the practitioner sufficient time for loading the cement into the cavity & sculpting it.

**Please read IFU for mixing, working and setting time details*



Impermeable Glass Bottle

Powder is packed in a pharma grade Glass bottle with tamper evident cap



Accurate Dispensing

Dropper with metered dose enables accurate dispensing with ideal powder to Liquid ratio using 1 scoop of powder & 1 drop of liquid



High Radiopacity

High Radio opacity allows for better Post-Operative Checks. This is due to the presence of Strontium in the Glass



TECHNICAL DATA

| | | |
|--------------|---------|----------|
| Mixing Time | @ 23° C | 0:45 sec |
| Working Time | @ 23° C | 2:00 min |
| Setting Time | @ 23° C | 3:30 min |

Packaging :

- 1 x 12 g Restore Plus powder bottle
- 1 x 7 g Restore Plus Liquid dropper
- Mixing pad
- Spoon

(Shelf life: 2 years from date of manufacture)

D- Tech Dental Technologies
62, Reality Warehousing
Wagholi, Pune-412207 INDIA
Url: www.dtechasia.com
ig: @dtech3d

For orders contact :
+91 9370145806
+91 9370145807
info@dtechasia.com

Proudly Made in India
for the World



CE MDD Certified
EN ISO13485:2016
ICMED Issue 2
REG NO.MFG/MD/2023/000135