

EMBRACE LC

Light Cure orthodontic adhesive bonding system

D Tech Product Newsletter

CE 0123

Optimal Bonding Strength

Embrace LC is a highly filled material that delivers optimum bond strength in all situations. Small particle filler for extra strength and yet easy to debond.

Tackiness

Embrace LC has the right level of tackiness. Ideal level of this property ensures that the cement holds the bracket in position and does not allow the bracket to drift during the adjustment of bracket position.


No solvent

Our unique formulation of primer is solvent free so there is no need to air dry the layer of primer before curing. This makes the priming less technique sensitive.

Perfect Debonding

Small particle filler for extra strength and yet easy to debond. It is essential that a bracket must debond in a manner that is safe for the enamel and at the same time not leaving behind too much cement for clean up.

d-tech



Translucency of paste

D Tech product is ideal to predicate M in translucency. This translucency indicates a perfect Refractive Index match of the fillers and resins which leads to a better curing throughout the cement. This is particularly important in bracket bonding because light exposure is only in the periphery of the bracket.

Packaging : 1x5g Adhesive Paste syringe
1x3ml Primer bottle
1x3g Ortho Etch syringe

D Tech Dental Technologies
Plot No 118/119, Sub Plot no 32
Lane No 7, Ramtekdi Industrial Estate
Hadapsar, Pune India 411028
Customer care- +91 9370145806
E mail: info@dtechasia.com
Url :www.dtechasia.com



CE MDD Certified
EN ISO13485:2016
ICMED Issue 2

Made in India
MFG LIC NO: MFG/MD/2024/000732



For orders contact :

+91 9370145806 +91 9370145807
+91 9370145808



info@dtechasia.com
ig: @dtech3d