

RESTORE +

High strength posterior Restorative cement

D Tech Product Newsletter

CE 0123

d-tech

Condensable

Premium quality ,condensable high strength posterior restorative cement for long term restorations with high compressive strength .



Comfort Set

Cement remains in condensable state up to 2 mins from start of mixing, giving the practitioner sufficient time for loading the cement into the cavity & sculpting it.

**Please read IFU for mixing, working and setting time details*



Impermeable Glass Bottle

Powder is packed in a pharma grade Glass bottle with tamper evident cap



Accurate Dispensing

Dropper with metered dose enables accurate dispensing with ideal powder to Liquid ratio using 1 scoop of powder & 1 drop of liquid



High Radiopacity

High Radio opacity allows for better Post-Operative Checks. This is due to the presence of Strontium in the Glass



TECHNICAL DATA

Mixing Time	@ 23° C	0:45 sec
Working Time	@ 23° C	2:00 min
Setting Time	@ 23° C	3:30 min

Packaging :

1 x 12 g Restore Plus powder bottle
 1 x 7 g Restore Plus Liquid dropper
 Mixing pad
 Spoon

(Shelf life: 2 years from date of manufacture)

D Tech Dental Technologies
 Plot No 118/119, Sub Plot no 32
 Lane No 7, Ramtekdi Industrial Estate
 Hadapsar , Pune India 411028
 Customer care- +91 9370145806
 E mail: info@dtechasia.com
 Url :www.dtechasia.com



CE MDD Certified
 EN ISO13485:2016
 ICMED Issue 2

Made in India
 MFG LIC NO: MFG/MD/2024/000732



For orders contact :



+91 9370145806 +91 9370145807
 +91 9370145808

info@dtechasia.com
 ig: @dtech3d